

CERTIFICATE OF ANALYSIS

Bovine Respiratory Syncytial Virus (BRSV)


FA Substrate Slide

Catalog No.:	SLD-IFA-BRSV
Lot:	P240225
Expiration:	08 April 2032

Quality Control Method: Indirect FA using BRSV positive control (catalog no. PC-IFA-BRSV), negative control (catalog no. NC-IFA-BRSV), and anti-bovine IgG_{1,2} affinity purified FITC conjugate (catalog no. CJ-F-BOVG-AP-10ML).

Specific Reaction: 3–4+ signal with the positive control and negative with the negative control. No background. There are 0–25 infected cells per high-power field.

Other Comments: NA

QC Manager: 
David Hyatt

Pattern Of Fluorescence: Granular cytoplasmic fluorescence with some small syncytia.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

11 April 2024