

CERTIFICATE OF ANALYSIS

Bluetongue Virus (BTV)

FA Substrate Slide

| Catalog No.: | SLD-IFA-BTV |
|--------------|---------------|
| Lot: | P210305-001 |
| Expiration: | 05 April 2025 |

Quality Control Method: Indirect FA using BTV positive control (catalog no. PC-IFA-BTV), negative

control (catalog no. NC-IFA-BTV), and anti-bovine IgG_{1,2} affinity purified

FITC conjugate (catalog no. CJ-F-BOVG-AP-10ML).

Specific Reaction: 2–4+ signal with the positive control, negative with the negative control,

and no background. There are 0-25 of infected cells per high-power field.

Other Comments: NA

QC Manager: 🔾

David Hyatt

Pattern Of Fluorescence: Multiple buckshot inclusions with dusty cytoplasmic fluorescence.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.