

CERTIFICATE OF ANALYSIS

Clostridium spp. 4-way

FA Substrate Slide

Catalog No.:	SLD-IFA-C4
Lot:	P180716-001
Expiration:	19 July 2022

Quality Control Method: Direct FA using conjugates: C. chauvoei (catalog no. CJ-F-CCO-1ML or

10ML), *C. novyi* (catalog no. CJ-F-CNO-1ML or 10ML), *C. septicum* (catalog no. CJ-F-CSE-1ML or 10ML), and *C. sordellii* (catalog no. CJ-F-CSO-1ML or

10ML).

Specific Reaction: C. chauvoei had 2-3+ signal, C. novyi had 1-3+ signal, C. septicum had 3-

4+ signal and *C. sordellii* had 2-3+ signal. There are greater than 10

organisms per high powered field.

Other Comments: NA

QC Manager: Chris Cay

QC Release Date: 20 July 2018

Pattern Of Fluorescence: Individual rod-shaped bacteria. Some clumped bacteria. Flagella

fluorescence alone is not indicative of a positive signal.

Technical Data Sheet Version: Version # 2