

CERTIFICATE OF ANALYSIS

Brucella canis

FA Substrate Slide

Catalog No.:	SLD-IFA-CB
Lot:	P160915-003
Expiration:	17 November 2020

Quality Control Method: IFA using *Brucella canis* positive control (catalog no. PC-IFA-CB), negative control (catalog no. NC-IFA-CB), and anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML).

Specific Reaction: 3+ signal with the positive control, negative with the negative control, no background noted with either control. Organisms are too numerous to count per high powered field.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 18 November 2016

Pattern Of Fluorescence: Brightly staining short bacilli against uniform background of ovine red blood cells.

Technical Data Sheet Version: Version #2

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\8-, 10-, 12-well Slides\Brucella canis\SLD-IFA-CB_P160915-003_201117.docx

18 November 2016