

CERTIFICATE OF ANALYSIS

Brucella canis

FA Substrate Slide

Catalog No.:	SLD-IFA-CB
Lot:	P190102-003
Expiration:	22 March 2023

Quality Control Method: IFA using *Brucella canis* positive control (catalog no. PC-IFA-CB), negative control (catalog no. NC-IFA-CB), and anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML).

Specific Reaction: 3+ signal with the positive control at neat, no background. Organisms are 50 to too numerous to count per high powered field.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 22 March 2019

Pattern of Fluorescence: Brightly staining short bacilli against uniform background of ovine red blood cells.

Technical Data Sheet Version: Version #2

FOR IN VITRO LABORATORY USE ONLY.

27 March 2019