

CERTIFICATE OF ANALYSIS

*Brucella canis*

FA Substrate Slide

Catalog No.:	SLD-IFA-CB
Lot:	P231221
Expiration:	26 December 2031

**Quality Control Method:** Indirect FA using *Brucella canis* positive control (catalog no. PC-IFA-CB), negative control (catalog no. NC-IFA-CB), and anti-canine IgG affinity purified FITC conjugate (catalog no. CJ-F-CANG-AP-10ML).

**Specific Reaction:** 3–4+ signal with the positive control and negative with the negative control. No background. Organisms per high-power field are too numerous to count.

**Other Comments:** NA

**QC Manager:**



David Hyatt

**Pattern Of Fluorescence:** Brightly staining short bacilli against uniform background of ovine red blood cells.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

03 January 2024