

CERTIFICATE OF ANALYSIS

Clostridium chauvoei

FA Substrate Slide

Catalog No.:	SLD-IFA-CCO
Lot:	P180430-002
Expiration:	04 June 2022

Quality Control Method: Direct FA using conjugates: *C. chauvoei* (catalog no. CJ-F-CCO-1ML or 10ML), *C. novyi* (catalog no. CJ-F-CNO-1ML or 10ML), *C. septicum* (catalog no. CJ-F-CSE-1ML or 10ML), and *C. sordellii* (catalog no. CJ-F-CSO-1ML or 10ML).

Specific Reaction: 3 – 4+ signal with *C. chauvoei*. Negative with trace background signal with *C. novyi*, *C. septicum*, and *C. sordellii*. Organisms are too numerous to count per high power field.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 04 June 2018

Pattern Of Fluorescence: Individual rod-shaped bacteria. Some clumped bacteria. Flagella fluorescence alone is not indicative of a positive signal.

Technical Data Sheet Version: Version # 2

FOR *IN VITRO* LABORATORY USE ONLY.

20 June 2018