

CERTIFICATE OF ANALYSIS

Clostridium chauvoei

FA Substrate Slide

Catalog No.:	SLD-IFA-CCO
Lot:	P180430-002
Expiration:	24 May 2026

Quality Control Method: Direct FA using conjugates: *C. chauvoei* (catalog no. CJ-F-CCO-1ML or 10ML), *C. novyi* (catalog no. CJ-F-CNO-1ML or 10ML), *C. septicum* (catalog no. CJ-F-CSE-1ML or 10ML), and *C. sordellii* (catalog no. CJ-F-CSO-1ML or 10ML).

Specific Reaction: 4+ signal with *C. chauvoei* conjugate. Trace–1+ signal with *C. novyi*, *C. septicum*, and *C. sordellii* conjugates. Organisms per high-power field are too numerous to count.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Individual rod-shaped bacteria. Some clumped bacteria. Flagella fluorescence alone is not indicative of a positive signal.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

26 May 2022