

CERTIFICATE OF ANALYSIS

Canine Coronavirus (CCV)

FA Control Slide

Catalog No.:	SLD-FAC-CCV
Lot:	P170508-001
Expiration:	21 May 2025

Quality Control Method: Direct FA using CCV conjugate (catalog no CJ-F-CCV-1ML or 10ML).

Specific Reaction: 2-4+ signal with trace-1+ background on positive well and negative with trace-1+ background on the negative well. There are 0-25 infected cells and/or syncytia per high powered field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern of Fluorescence: Large fluorescent syncytia with some single cells in between.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\2-well Slides\CCV\SLD-IFA-CCV_P170508-001_250521.docx

02 June 2021