

CERTIFICATE OF ANALYSIS

Canine Coronavirus (CCV)

FA Substrate Slide

Catalog No.:	SLD-IFA-CCV
Lot:	P200520-003
Expiration:	04 June 2024

Quality Control Method: Indirect FA using CCV antiserum (catalog no PAB-CCV) and anti-porcine affinity purified IgG FITC conjugate (catalog no CJ-F-PORG-AP-10ML).

Specific Reaction: 2-3+ signal with the antiserum at 1/400 and negative with the diluent control with no background. There are 0-50 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Large fluorescent syncytia with some single cells in between.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.

10 June 2020