

CERTIFICATE OF ANALYSIS

## Canine Coronavirus (CCV)

FA Substrate Slide

Catalog No.:	SLD-IFA-CCV
Lot:	P220506-001
Expiration:	10 June 2026

**Quality Control Method:** Indirect FA using FIP-1 positive control (catalog no. PC-IFA-FIP1), negative control (catalog no NC-IFA-FIP1), and anti-feline IgG affinity purified FITC conjugate (catalog no. CJ-F-FELG-AP-10ML).

**Specific Reaction:** 2–4+ signal with the positive control and negative with the negative control. No background. There are 0–50 infected cells per high-power field.

**Other Comments:** NA

**QC Manager:**



David Hyatt

**Pattern Of Fluorescence:** Large fluorescent syncytia with some single cells in between.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.

14 June 2022