

**CERTIFICATE OF ANALYSIS**

# Canine Distemper Virus (CDV)

FA Substrate Slide

Catalog No.:	SLD-IFA-CDV
Lot:	P180615-002
Expiration:	13 August 2022

**Quality Control Method:** Indirect FA using CDV IgG positive control (catalog no. PC-IFA-CDV-G), anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML), CDV IgM positive control (catalog no. PC-IFA-CDV-M), and anti-canine IgM affinity purified FITC conjugate (catalog no. CJ-F-CANM-AP-1ML or 10ML).

**Specific Reaction:** 2-4+ signal with the IgG positive control. 1-2+ with the IgM positive control. Negative with the diluent and no background. There are 0-50 infected cells per high powered field.

**Other Comments:** NA

**QC Manager:** Chris Cay



**QC Release Date:** 20 August 2018

**Pattern Of Fluorescence:** Individual cells with inclusion bodies and some plasma membrane fluorescence.

**Technical Data Sheet Version:** Version # 2

FOR *IN VITRO* LABORATORY USE ONLY.

27 August 2018