

CERTIFICATE OF ANALYSIS

Canine Distemper Virus (CDV)

FA Substrate Slide

Catalog No.:	SLD-IFA-CDV
Lot:	P180615-003
Expiration:	28 August 2022

Quality Control Method: Indirect FA using CDV IgG positive control (catalog no. PC-IFA-CDV-G), anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML), CDV IgM positive control (catalog no. PC-IFA-CDV-M), and anti-canine IgM affinity purified FITC conjugate (catalog no. CJ-F-CANM-AP-1ML or 10ML).

Specific Reaction: 3 – 4+ signal with IgG positive control. 2 – 4+ signal with IgM positive control. Negative with diluent and no background. There are 0 – 50 infected cells per high power field.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 31 August 2018

Pattern Of Fluorescence: Individual cells with inclusion bodies and some plasma membrane fluorescence.

Technical Data Sheet Version: Version # 2

FOR *IN VITRO* LABORATORY USE ONLY.

17 October 2018