

CERTIFICATE OF ANALYSIS

Canine Distemper Virus (CDV)

FA Substrate Slide

| | |
|--------------|-----------------|
| Catalog No.: | SLD-IFA-CDV |
| Lot: | P211005-001 |
| Expiration: | 28 January 2026 |

Quality Control Method: Indirect FA using CDV IgG positive control (catalog no. PC-IFA-CDV-G), anti-canine IgG affinity purified FITC conjugate (catalog no. CJ-F-CANG-AP-10ML), CDV IgM positive control (catalog no. PC-IFA-CDV-M), and anti-canine IgM affinity purified FITC conjugate (catalog no. CJ-F-CANM-AP-10ML).

Specific Reaction: 3–4+ signal, trace–1+ background with the IgG positive control. 2–4+ signal, trace–1+ background with the IgM positive control. Negative, trace background with the diluent control. There are 0–50 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Individual cells with inclusion bodies and some plasma membrane fluorescence.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

08 February 2022