

CERTIFICATE OF ANALYSIS

Canine Distemper Virus (CDV)

FA Substrate Slide

Catalog No.:	SLD-IFA-CDV
Lot:	P211005-002
Expiration:	09 February 2026

Quality Control Method: Indirect FA using CDV IgG positive control (catalog no. PC-IFA-CDV-G), anti-canine IgG affinity purified FITC conjugate (catalog no. CJ-F-CANG-AP-10ML), CDV IgM positive control (catalog no. PC-IFA-CDV-M), and anti-canine IgM affinity purified FITC conjugate (catalog no. CJ-F-CANM-AP-10ML).

Specific Reaction: 4+ signal, trace background with the IgG positive control. 2-4+ signal, trace background with the IgM positive control. Negative with the diluent control. There are 0-50 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Individual cells with inclusion bodies and some plasma membrane fluorescence.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

14 February 2022