

CERTIFICATE OF ANALYSIS

Canine Herpesvirus Type 1 (CHV-1)

FA Substrate Slide

Catalog No.:	SLD-IFA-CHV
Lot:	P160303-002
Expiration:	28 January 2024

Quality Control Method: Indirect FA using CHV-1 positive control (catalog no. PC-IFA-CHV) and anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML).

Specific Reaction: 2-4+ signal with trace - 1 background with the positive control and negative with the diluent control. There are 0-20 of infected cells or plaques per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Cytoplasmic and nuclear fluorescence in rounded cells with occasional small syncytia.

Technical Data Sheet Version: Version # 1

FOR IN VITRO LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\8-, 10-, 12-well Slides\CHV\SLD-IFA-CHV_P160303-002_240128.docx

3 February 2020