

## **CERTIFICATE OF ANALYSIS**

## Canine Herpesvirus Type 1 (CHV-1)

FA Substrate Slide

Catalog No.:	SLD-IFA-CHV
Lot:	P200713-001
Expiration:	31 July 2024

Quality Control Method: Indirect FA using CHV-1 positive control (catalog no. PC-IFA-CHV) and

anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-10ML)

Specific Reaction: 2–4+ signal with the positive control and negative with the diluent control

with trace-1+ background. Infected cells per high-power field are too

numerous to count.

Other Comments: NA

QC Manager:

David Hyatt

Pattern Of Fluorescence: Cytoplasmic and nuclear fluorescence in rounded cells with occasional

small syncytia

Technical Data Sheet Version: 1