

CERTIFICATE OF ANALYSIS

## Canine Herpesvirus Type 1 (CHV-1)

FA Substrate Slide

Catalog No.:	SLD-IFA-CHV
Lot:	P240119
Expiration:	26 February 2032

**Quality Control Method:** Indirect FA using CHV-1 positive control (catalog no. PC-IFA-CHV) and anti-canine IgG affinity purified FITC conjugate (catalog no. CJ-F-CANG-AP-10ML).

**Specific Reaction:** 3–4+ signal, trace–1+ background with the positive control and negative with the negative control. Infected cells per high-power field are too numerous to count.

**Other Comments:** NA

**QC Manager:**



David Hyatt

**Pattern Of Fluorescence:** Cytoplasmic and nuclear fluorescence in rounded cells with occasional small syncytia.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

15 March 2024