

CERTIFICATE OF ANALYSIS

Canine Parvovirus (CPV)

FA Substrate Slide

Catalog No.:	SLD-IFA-CPV
Lot:	P181203-002
Expiration:	18 February 2023

Quality Control Method: Indirect FA using CPV IgG positive control (catalog no. PC-IFA-CPV-G),

CPV IgM positive control (catalog no. PC-IFA-CPV-M), CPV In-house negative control, anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML), and anti-canine IgM FITC affinity purified conjugate (catalog

no. CJ-F-CANM-AP-1ML or 10ML).

Specific Reaction: 2-4+ signal with the IgG positive control, no background. 1-3+ signal with

the IgM positive control, no background. 0-50 infected cells per high

powered field.

Other Comments: NA

QC Manager: Chris Cay

QC Release Date: 19 February 2019

Pattern of Fluorescence: Single cells with cytoplasmic and nuclear fluorescence.

Technical Data Sheet Version: Version #2