

CERTIFICATE OF ANALYSIS

Ehrlichia canis

FA Substrate Slide

Catalog No.:	SLD-IFA-EC
Lot:	P180803-003
Expiration:	11 October 2022

Quality Control Method: Indirect FA using *Ehrlichia canis* positive control (catalog no. PC-IFA-EC), negative control (catalog no. NC-IFA-EC), and anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML).

Specific Reaction: 3-4+ signal with the positive control with no background and negative with the negative control with no background. There are 0 to 50 infected cells per high-power field.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 15 October 2018

Pattern Of Fluorescence: Brightly fluorescent cytoplasmic inclusion bodies (morulae), clustered elementary bodies (initial bodies), and free elementary bodies.

Technical Data Sheet Version: Version # 2

FOR IN VITRO LABORATORY USE ONLY.

17 October 2018