

CERTIFICATE OF ANALYSIS

*Ehrlichia canis*

FA Substrate Slide

Catalog No.:	SLD-IFA-EC
Lot:	P181003-001
Expiration:	31 October 2022

Quality Control Method: Indirect FA using *Ehrlichia canis* positive control (catalog no. PC-IFA-EC), negative control (catalog no. NC-IFA-EC), and anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML).

Specific Reaction: 2-4+ signal with the positive control with trace background and negative with the negative control with trace background. There is an average of 0 to 50 infected cells per high power field.

Other Comments: NA

QC Manager: Chris Cay



QC Release Date: 01 November 2018

Pattern Of Fluorescence: Brightly fluorescent cytoplasmic inclusion bodies (morulae), clustered elementary bodies (initial bodies), and free elementary bodies.

Technical Data Sheet Version: Version # 2

FOR IN VITRO LABORATORY USE ONLY.

15 November 2018