

CERTIFICATE OF ANALYSIS

Ehrlichia canis

FA Substrate Slide

| | |
|--------------|-------------|
| Catalog No.: | SLD-IFA-EC |
| Lot: | P230601-001 |
| Expiration: | 25 May 2031 |

Quality Control Method: Indirect FA using *Ehrlichia canis* positive control (catalog no. PC-IFA-EC), negative control (catalog no. NC-IFA-EC), and anti-canine IgG affinity purified FITC conjugate (catalog no. CJ-F-CANG-AP-10ML).

Specific Reaction: 2–4+ signal with the positive control and negative with the negative control. No background. There are more than 50 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Brightly fluorescent cytoplasmic inclusion bodies (morulae), clustered elementary bodies (initial bodies), and free elementary bodies.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.

02 June 2023