

CERTIFICATE OF ANALYSIS

Feline Calicivirus (FCV)

FA Substrate Slide

Catalog No.:	SLD-IFA-FCV
Lot:	P170127-001
Expiration:	30 March 2021

Quality Control Method: IFA using FCV positive control (catalog no. PC-IFA-FCV), FCV negative control (catalog no. NC-IFA-FCV), and anti-feline IgG affinity purified conjugate (catalog no. CJ-F-FELG-AP-1ML or 10ML).

Specific Reaction: 2-4+ signal with the positive control at neat, no background. Negative with the negative control, no background. There are 0-30 infected cells per high-power field.

Other Comments: NA

QC Manager:  Will Thurman

QC Release Date: 31 March 2017

Pattern of Fluorescence: Diffuse and granular cytoplasmic fluorescence with brightly fluorescing nuclear membranes, also brightly fluorescing individual nuclei.

Technical Data Sheet Version: Version #2

FOR IN VITRO LABORATORY USE ONLY.

H:\Quality\COA (VMRD)\8-, 10-, 12-well Slides\FCV\SLD-IFA-FCV_P170127-001_210330.docx

19 April 2017