

CERTIFICATE OF ANALYSIS

Feline Calicivirus (FCV)

FA Substrate Slide

Catalog No.:	SLD-IFA-FCV
Lot:	P220902-001
Expiration:	17 October 2030

Quality Control Method: Indirect FA using FCV positive control (catalog no. PC-IFA-FCV), FCV negative control (catalog no. NC-IFA-FCV), and anti-feline IgG affinity purified conjugate (catalog no. CJ-F-FELG-AP-10ML).

Specific Reaction: 2–4+ signal with the positive control and negative with the negative control. No background. There are 0–25 infected cells per high-power field.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Diffuse and granular cytoplasmic fluorescence with brightly fluorescing nuclear membranes, also brightly fluorescing individual nuclei.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.

25 October 2022