

**CERTIFICATE OF ANALYSIS**

# Feline Infectious Peritonitis Virus Type 1 (FIP-1)

FA Substrate Slide

Catalog No.:	SLD-IFA-FIP1
Lot:	P180627-001
Expiration:	26 September 2022

**Quality Control Method:** Indirect FA using FIP-1 positive control (catalog no. PC-IFA-FIP1), negative control (catalog no NC-IFA-FIP1), and anti-feline IgG FITC affinity purified conjugate (catalog no. CJ-F-FELG-AP-1ML or 10ML).

**Specific Reaction:** 2 – 4+ signal with the positive control and negative with the negative control, no background. 0 – 50 infected cells per high power field.

**Other Comments:** NA

**QC Manager:** Chris Cay



**QC Release Date:** 03 October 2018

**Pattern Of Fluorescence:** Multiple syncytia, plaques and some individual cells with granular cytoplasmic fluorescence.

**Technical Data Sheet Version:** Version # 2

FOR *IN VITRO* LABORATORY USE ONLY.

18 October 2018