

CERTIFICATE OF ANALYSIS

# Feline Infectious Peritonitis Virus Type 1 (FIP-1)


FA Substrate Slide

|              |               |
|--------------|---------------|
| Catalog No.: | SLD-IFA-FIP1  |
| Lot:         | P200215-001   |
| Expiration:  | 12 March 2024 |

**Quality Control Method:** Indirect FA using FIP-1 positive control (catalog no. PC-IFA-FIP1), negative control (catalog no NC-IFA-FIP1), and anti-feline IgG FITC affinity purified conjugate (catalog no. CJ-F-FELG-AP-10ML).

**Specific Reaction:** 3-4+ signal with the positive control and trace with the negative control with no background. There are 0-100 infected cells per high-power field.

**Other Comments:** NA

**QC Manager:**   
David Hyatt

**Pattern Of Fluorescence:** Multiple syncytia, plaques and some individual cells with granular cytoplasmic fluorescence.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

13 April 2020