

## **CERTIFICATE OF ANALYSIS**

## Feline Infectious Peritonitis Virus Type 1 (FIP-1)

FA Substrate Slide

Catalog No.:	SLD-IFA-FIP1
Lot:	P220913-002
Expiration:	07 December 2030

Quality Control Method: Indirect FA using FIP-1 positive control (catalog no. PC-IFA-FIP1), negative

control (catalog no NC-IFA-FIP1), and anti-feline IgG affinity purified FITC

conjugate (catalog no. CJ-F-FELG-AP-10ML).

**Specific Reaction**: 2–4+ signal with the positive control and negative with the negative control.

No background. There are 0–50 infected cells per high-power field.

Other Comments: NA

QC Manager: 🔾

David Hyatt

Pattern Of Fluorescence: Multiple syncytia, plaques and some individual cells with granular

cytoplasmic fluorescence.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.