

CERTIFICATE OF ANALYSIS

Feline Herpesvirus Type 1 (FHV-1/FVR)

FA Substrate Slide

Catalog No.:	SLD-IFA-FVR
Lot:	P141118-004
Expiration:	12 January 2019

Quality Control Method: Indirect FA using FHV-1/FVR positive control (catalog no. PC-IFA-FVR), and

anti-feline IgG FITC affinity purified conjugate (catalog no. CJ-F-FELG-AP-

1ML or 10ML).

Specific Reaction: 2-4+ fluorescence on the positive control and negative on the diluent

control with no background. There are 0-30 plagues or infected cells per

high-power field.

NA Other Comments:

> Joney Buson Torrey Beeson QC Manager:

QC Release Date: 13 January 2015

Pattern Of Fluorescence: Individual rounded cells with cytoplasmic and nuclear fluorescence. Some

small syncytia.

Technical Data Sheet Version: Version # 1

www.vmrd.com

PO Box 502 P: 509.334.5815 techserve@vmrd.com Pullman, WA 99163 USA F: 509.332.5356