

CERTIFICATE OF ANALYSIS

Feline Herpesvirus Type 1 (FHV-1/FVR)

FA Substrate Slide

| | |
|--------------|----------------|
| Catalog No.: | SLD-IFA-FVR |
| Lot: | P200629-001 |
| Expiration: | 03 August 2024 |

Quality Control Method: Indirect FA using FHV-1/FVR positive control (catalog no. PC-IFA-FVR), and anti-feline IgG FITC affinity purified conjugate (catalog no. CJ-F-FELG-AP-10ML).

Specific Reaction: 2–4+ signal with the positive control and negative with the diluent control with no background. There are 0–50 infected cells per high-power field.

Other Comments: NA

QC Manager: 
David Hyatt

Pattern Of Fluorescence: Individual rounded cells with cytoplasmic and nuclear fluorescence. Some small syncytia.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

10 August 2020