

## **CERTIFICATE OF ANALYSIS**

## Borrelia burgdorferi (Lyme Disease)

FA Substrate Slide

Catalog No.:	SLD-IFA-LD
Lot:	P171006-010
Expiration:	6 February 2022

Quality Control Method: Indirect FA using Borrelia burgdorferi positive control (catalog no. PC-IFA-

LD), negative control (catalog no. NC-IFA-LD), and anti-canine IgG FITC

conjugate (catalog no. CJ-F-CANG-1ML or 10ML).

Specific Reaction: 2 - 4+ signal with the positive control and negative with the negative

control, no background. There are 100 to too numerous to count

organisms per high-power field.

Other Comments: NA

QC Manager: Wile Thurman Will Thurman

QC Release Date: 12 February 2018

Pattern Of Fluorescence: Spirochetes with uniform surface fluorescence

Technical Data Sheet Version: Version # 2