

CERTIFICATE OF ANALYSIS

Borrelia burgdorferi (Lyme Disease)

FA Substrate Slide

Catalog No.:	SLD-IFA-LD
Lot:	P190909-003
Expiration:	05 March 2024

Quality Control Method: Indirect FA using Borrelia burgdorferi positive control (catalog no. PC-IFA-

LD), negative control (catalog no. NC-IFA-LD), and anti-canine IgG FITC

conjugate (catalog no. CJ-F-CANG-10ML).

Specific Reaction: 2-4+ signal with the positive control with trace background. Negative with

the negative control. There are 0-50 organisms per high-power field.

Other Comments: NA

David Hyatt

Pattern Of Fluorescence: Spirochetes with uniform surface fluorescence.

Technical Data Sheet Version: 2

FOR IN VITRO LABORATORY USE ONLY.