

CERTIFICATE OF ANALYSIS

## Porcine Circovirus Type 2 (PCV-2)

FA Substrate Slide

Catalog No.:	SLD-IFA-PCV2
Lot:	P180110-004
Expiration:	02 December 2025

**Quality Control Method:** Indirect FA using PCV-2 positive control (catalog no. PC-IFA-PCV2), negative control (catalog no. NC-IFA-PCV2), and anti-porcine IgG affinity purified conjugate (catalog no. CJ-F-PORG-AP-1ML or 10ML).

**Specific Reaction:** 3–4+ signal with the positive control and negative with the negative control with trace background. There are 0–25 infected cells per high-power field.

**Other Comments:** NA

**QC Manager:**



David Hyatt

**Pattern Of Fluorescence:** Intense whole nuclear and variable sized but discrete cytoplasmic inclusions

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Master Documents\COA (VMRD)\8-, 10-, 12-well Slides\PCV2\SLD-IFA-PCV2\_P180110-004\_251202.docx

03 December 2021