

CERTIFICATE OF ANALYSIS

Porcine Epidemic Diarrhea Virus (PEDV)


FA Substrate Slide

Catalog No.:	SLD-IFA-PEDV
Lot:	P150501-003
Expiration:	23 March 2031

Quality Control Method: Indirect FA using PEDV positive control (catalog no. PC-IFA-PEDV), negative control (catalog no. NC-IFA-PEDV), and anti-porcine IgG affinity purified conjugate (catalog no. CJ-F-PORG-AP-10ML).

Specific Reaction: 1–3+ signal with the positive control and negative with the negative control. No background. There are 0–25 infected cells per high-power field.

Other Comments: NA

QC Manager: 

David Hyatt

Pattern Of Fluorescence: Brightly fluorescing irregular cytoplasmic inclusions, mainly in single cells.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

4 April 2023