

CERTIFICATE OF ANALYSIS

Trypanosoma cruzi

FA Substrate Slide

Catalog No.:	SLD-IFA-TC
Lot:	P220705-001
Expiration:	02 December 2030

Quality Control Method: Indirect FA using *Trypanosoma cruzi* antibody positive human serum and anti-human IgG FITC conjugate.

Specific Reaction: 3–4+ signal with positive serum and negative with the diluent control. No background. Organisms per high-power field are too numerous to count.

Other Comments: NA

QC Manager:



David Hyatt

Pattern Of Fluorescence: Brightly fluorescing elongated trypanosome organisms with clearly defined flagella.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

07 December 2022