

CERTIFICATE OF ANALYSIS

*Trypanosoma cruzi*

FA Substrate Slide

Catalog No.:	SLD-IFA-TC
Lot:	P231025
Expiration:	22 November 2031

**Quality Control Method:** Indirect FA using *Trypanosoma cruzi* antibody positive human serum and anti-human IgG FITC conjugate.

**Specific Reaction:** 3–4+ signal with positive serum and negative with the diluent control. No background. Organisms per high-power field are too numerous to count.

**Other Comments:** NA

**QC Manager:**



David Hyatt

**Pattern Of Fluorescence:** Brightly fluorescing elongated trypanosome organisms with clearly defined flagella.

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

05 December 2023