

CERTIFICATE OF ANALYSIS

Transmissible Gastroenteritis Virus (TGEV)

FA Control Slide

Catalog No.:	SLD-FAC-TGE
Lot:	P170207-001
Expiration:	17 December 2024

Quality Control Method: Direct FA using VMRD, Inc. TGE Direct FA Conjugate (catalog no. CJ-F-

TGE-10ML).

Specific Reaction: 3–4+ signal with trace background on the positive well and negative with

trace background on the negative well. There are 0–10 infected cells

and/or syncytia per high-power field.

Other Comments: NA

QC Manager: (

David Hyatt

Pattern of Fluorescence: Moderate to large syncytia with granular and smooth cytoplasmic

P: 509.334.5815

F: 509.332.5356

fluorescence.

Technical Data Sheet Version: 2

techserve@vmrd.com

www.vmrd.com