

CERTIFICATE OF ANALYSIS

Toxoplasma gondii

FA Substrate Slide

Catalog No.:	SLD-IFA-TOXO
Lot:	P190228-001
Expiration:	17 April 2023

Quality Control Method: Indirect FA using *Toxoplasma gondii* IgG positive control (catalog no. PC-IFA-TOXO-FEL-G), negative control (catalog no. NC-IFA-TOXO-FEL), anti-feline IgG FITC affinity purified conjugate (catalog no. CJ-F-FELG-AP-1ML or 10ML).

Specific Reaction: 3+ signal with the IgG positive control at neat, no background. Negative with the negative control with 1-2+ apical staining. Organisms were 50 to too numerous to count per high powered field.

Other Comments: When using this slide with canine serum it is recommended to use VMRD's anti-canine affinity purified conjugate (catalog no. CJ-F-CANG-AP-1ML or 10ML) and FA Serum Diluting Buffer with 1% BSA pH 7.2 (catalog no. FASDB-100ML) as a diluent.

QC Manager: Chris Cay



QC Release Date: 18 April 2019

Pattern of Fluorescence: Organisms with bright diffuse or peripheral fluorescence.

Technical Data Sheet Version: Version # 2

FOR IN VITRO LABORATORY USE ONLY.

22 April 2019