

CERTIFICATE OF ANALYSIS

*Toxoplasma gondii*


FA Substrate Slide

Catalog No.:	SLD-IFA-TOXO
Lot:	P220415-001
Expiration:	01 June 2026

**Quality Control Method:** Indirect FA using *Toxoplasma gondii* IgG positive control (catalog no. PC-IFA-TOXO-FEL-G), negative control (catalog no. NC-IFA-TOXO-FEL), anti-feline IgG affinity purified FITC conjugate (catalog no. CJ-F-FELG-AP-10ML).

**Specific Reaction:** 3–4+ signal, no background with the positive control. Negative with 1–2+ apical staining with the negative control. Organisms per high-power field were too numerous to count.

**Other Comments:** When using *Toxoplasma gondii* systems it is recommended to use FA Serum Diluting Buffer with 1% BSA pH 7.2 (catalog no. FASDB-100ML) as a diluent.

**QC Manager:**   
David Hyatt

**Pattern Of Fluorescence:** Organisms with bright diffuse or peripheral fluorescence.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

07 June 2022