

CERTIFICATE OF ANALYSIS

*Toxoplasma gondii*


FA Substrate Slide

Catalog No.:	SLD-IFA-TOXO
Lot:	P220415-004
Expiration:	08 July 2030

**Quality Control Method:** Indirect FA using *Toxoplasma gondii* IgG positive control (catalog no. PC-IFA-TOXO-FEL-G), negative control (catalog no. NC-IFA-TOXO-FEL), and anti-feline IgG affinity purified FITC conjugate (catalog no. CJ-F-FELG-AP-10ML).

**Specific Reaction:** 2–4+ signal, no background with the positive control. Negative with trace apical staining with the negative control. Organisms per high-power field are too numerous to count

**Other Comments:** When using this slide with canine serum it is recommended to use VMRD's anti-canine affinity purified conjugate (catalog no. CJ-F-CANG-AP-10ML) and FA Serum Diluting Buffer with 1% BSA pH 7.2 (catalog no. FASDB-100ML) as a diluent.

QC Manager: 

David Hyatt

**Pattern Of Fluorescence:** Organisms with bright diffuse or peripheral fluorescence

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

21 July 2022