

**CERTIFICATE OF ANALYSIS**

## Bluetongue Virus (BTV)

FA Substrate Slide

Catalog No.:	SLD-IFA-BTV
Lot:	P141107-004
Expiration:	20 November 2018

**Quality Control Method:** Indirect FA using BTV positive control (catalog no. PC-IFA-BTV), negative control (catalog no. NC-IFA-BTV), and anti-bovine IgG<sub>1,2</sub> affinity purified FITC conjugate (catalog no. CJ-F-BOVG-AP-1ML or 10ML).

**Specific Reaction:** 2-4+ fluorescence on the positive control and negative on the negative control with no background. There are 10 to too numerous to count infected cells per high-power field.

**Other Comments:** This lot of slides has considerably higher infection rate than previous lots.

**QC Manager:** *Torrey Beeson* Torrey Beeson

**QC Release Date:** 21 November 2014

**Pattern Of Fluorescence:** Multiple buckshot inclusions with dusty cytoplasmic fluorescence.

**Technical Data Sheet Version:** Version # 1

**FOR IN VITRO LABORATORY USE ONLY.**

H:\Quality VMRD\COAs and TDS\COA (VMRD)\8-, 10-, 12-well Slides\BTV\SLD-IFA-BTV\_P141107-004\_182011.docx

02 December 2014