

CERTIFICATE OF ANALYSIS

Canine Distemper Virus (CDV)

FA Substrate Slide

Catalog No.:	SLD-IFA-CDV
Lot:	P141118-005
Expiration:	16 December 2018

Quality Control Method: Indirect FA using CDV IgG positive control (catalog no. PC-IFA-CDV-G), anti-canine IgG FITC conjugate (catalog no. CJ-F-CANG-1ML or 10ML), CDV IgM positive control (catalog no. PC-IFA-CDV-M), and anti-canine IgM affinity purified FITC conjugate (catalog no. CJ-F-CANM-AP-1ML or 10ML).

Specific Reaction: 3-4+ fluorescence with the IgG positive control, 1-3+ with the IgM positive control and negative with the diluent control, no background. There are 0-50 infected cells per high-power field.

Other Comments: NA

QC Manager: *Torrey Beeson* Torrey Beeson

QC Release Date: 17 December 2014

Pattern Of Fluorescence: Individual cells with inclusion bodies and some plasma membrane fluorescence.

Technical Data Sheet Version: Version # 1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality VMRD\COAs and TDS\COA (VMRD)\8-, 10-, 12-well Slides\CDV\SLD-IFA-CDV_P141118-005_181612.docx

07 January 2015