

CERTIFICATE OF ANALYSIS

Rabies Recombinant Nucleoprotein (rNP)

FA Control Slide

Catalog No.:	SLD-FAC-RAB
Lot:	P141107-001
Expiration:	14 November 2018

Quality Control Method: Direct FA using Rabies FITC in-house conjugate.

Specific Reaction: 2-4+ fluorescence on the positive well and negative on the negative well

with no background. There are 0-20 infected cells per high-power field.

Other Comments: NA

QC Manager: Joney Busan Torrey Beeson

QC Release Date: 18 November 2014

Pattern Of Fluorescence: Smooth and particulate cytoplasmic fluorescence of single and

multinucleated cells.

Technical Data Sheet Version: Version # 1

P: 509.334.5815 F: 509.332.5356 techserve@vmrd.com www.vmrd.com