



VMRD, Inc.

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: vmrd@vmrd.com

Web site: <http://www.vmrd.com>

Certificate of Analysis


EQUINE IgG
AFFINITY PURIFIED
FITC ANTI-IMMUNOGLOBULIN
CONJUGATE

CATALOG NO.: 043-AP-10

VOLUME: 10 ml

LOT: P030122002-020605

EXPIRATION: 6 February 2005



DESCRIPTION: Anti-Equine IgG polyclonal antiserum, affinity purified, conjugated to fluorescein isothiocyanate. Caprine origin. Ready for use.

QUALITY CONTROL METHOD: IFA on VMRD, Inc. Equine Rhinopneumonitis Virus/Equine Herpesvirus Type I (ERV/EHV-1) Substrate Slides (catalog no. 210-88-10-ERV), Equine Origin ERV/EHV-1 Positive Control (catalog no. 211-P-ERV-EQ).

Specific Reaction: 4+ on positive cells; negative on negative cells.

Other Reactions or Comments: Trace background.

STORAGE: This conjugate is provided in liquid form and should be stored at 4-8°C. DO NOT FREEZE! If conjugate becomes cloudy it should be discarded. This conjugate contains 0.09% sodium azide as a preservative.

FOR IN VITRO LABORATORY USE ONLY.

WARRANTY: VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.