



VMRD, Inc.

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: vtech@vmrd.com

Web site: <http://www.vmrd.com>

Certificate of Analysis


FELINE IgM **FITC ANTI-IMMUNOGLOBULIN** **CONJUGATE**

CATALOG NO.: 056-10

VOLUME: 10 ml

LOT: V01071-082304

EXPIRATION: 23 August 2004



DESCRIPTION: Anti-Feline IgM polyclonal antiserum conjugated to fluorescein isothiocyanate. Heavy chain specific. Caprine origin. Affinity purified. Ready for use.

QUALITY CONTROL METHOD: Indirect FA using VMRD *Toxoplasma gondii* 12-well substrate slide (catalog no. 210-88-12-TOXO) to detect binding of feline anti-*Toxoplasma* IgM; *Toxoplasma* IgM positive and negative controls (internal control samples).

Specific Reaction: 4+ positive fluorescence on the whole cells. Negative control shows 4+ fluorescence only in apical region of the *Toxoplasma* cells.

Other Reactions or Comments: No background.

STORAGE: This conjugate is provided in liquid form and should be stored at 4-8°C. **DO NOT FREEZE!** If conjugate becomes cloudy it should be discarded. This conjugate contains 0.09% sodium azide as a preservative.

FOR RESEARCH AND INVESTIGATIONAL USE ONLY.

WARRANTY: VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.