



**VMRD, Inc.**

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: [vtech@vmrd.com](mailto:vtech@vmrd.com)

Web site: <http://www.vmrd.com>

## Certificate of Analysis

### **GOAT IgG**


**FITC ANTI-IMMUNOGLOBULIN  
CONJUGATE**

**CATALOG NO.:** 075-10

**VOLUME:** 10 ml

**LOT:** V02015-051306

**EXPIRATION DATE:** 13 May 2006



**DESCRIPTION:** Anti-Goat IgG polyclonal antiserum conjugated to fluorescein isothiocyanate. Rabbit origin. Ready for use.

**QUALITY CONTROL METHOD:** Indirect FA using VMRD Equine Rhinopneumonitis Virus/Equine Herpesvirus Type 1 (ERV/EHV-1) Substrate Slides (catalog no. 210-88-10-ERV), Goat Anti-ERV Antiserum (catalog no. 210-70-ERV); Bovine Viral Diarrhea Virus (BVDV) Substrate Slides (catalog no. 210-88-10-BVD), Goat Anti-BVDV Antiserum (catalog no. 210-70-BVD).

**Specific Reaction:** >4+ positive on ERV cells; 4+ positive on BVDV cells; negative on negative cells.

**Other Reactions or Comments:** Conjugate should be centrifuged before use to remove excess artifacts.

**STORAGE:** This conjugate is provided in liquid form and should be stored at 4-8°C. DO NOT FREEZE! If conjugate becomes cloudy it should be discarded. This conjugate contains 0.09% sodium azide as a preservative.

FOR RESEARCH AND INVESTIGATIONAL USE ONLY.

**WARRANTY:** VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.