



**VMRD**

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: [vmrd@vmrd.com](mailto:vmrd@vmrd.com)

Web site: <http://www.vmrd.com>

## Certificate of Analysis

### **LLAMA IgG**


### **FITC ANTI-IMMUNOGLOBULIN CONJUGATE**

**CATALOG NO.:** 066-10

**VOLUME:** 10 ml

**LOT:** V02069

**EXPIRATION DATE:** 21 June 2008



**DESCRIPTION:** Anti-Llama IgG polyclonal antiserum conjugated to fluorescein isothiocyanate. Affinity purified. Reacts with heavy and light chains. Goat origin. Ready for use.

**QUALITY CONTROL METHOD:** Indirect FA using VMRD, Inc. Equine Rhinopneumonitis Virus (ERV) FA Substrate Slides (catalog no. 210-88-10-ERV), ERV Llama Positive Control (catalog no. 211-P-ERV-LL).

**Specific Reaction:** 2-3+ positive on ERV Positive Control; negative on ERV negative cells.

**Other Reactions or Comments:** 1+ background at neat.

**STORAGE:** This conjugate is provided in liquid form and should be stored at 4-8°C. **DO NOT FREEZE!** If conjugate becomes cloudy it should be discarded. This conjugate contains 0.09% sodium azide as a preservative.

FOR *IN VITRO* LABORATORY USE ONLY.

**WARRANTY:** VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.