



**VMRD**

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: [vmrd@vmrd.com](mailto:vmrd@vmrd.com)

Web site: <http://www.vmrd.com>

## Certificate of Analysis

### **MOUSE IgG** **FITC ANTI-IMMUNOGLOBULIN** **CONJUGATE**

**CATALOG NO.:** 095-10

**VOLUME:** 10 ml

**LOT:** P080208-008

**EXPIRATION:** 18 February 2010



**DESCRIPTION:** Anti-Mouse IgG polyclonal antiserum conjugated to fluorescein isothiocyanate. IgG specific by immunoelectrophoresis (IEP). Rabbit origin. Affinity purified. Ready for use.

**QUALITY CONTROL METHOD:** Indirect FA using VMRD Canine Distemper Virus (CDV) 10-well substrate slide (catalog no. SLD-IFA-CDV), mouse origin CDV-specific monoclonal antibody (catalog no. CDV-NP-0.1), mouse origin Canine Parainfluenza Virus-specific monoclonal antibody (catalog no. CPI-A-CA-0.1).

**Specific Reaction:** 4+ positive on CDV-NP diluted to 5 µg/ml in VMRD Conjugate Diluting Buffer (catalog no. 210-91-CB).

**Other Reactions or Comments:** Trace 1+ background.

**STORAGE:** This conjugate is provided in liquid form and should be stored at 4-8°C. DO NOT FREEZE! If conjugate becomes cloudy it should be discarded. This conjugate contains 10 ppm ProClin 300 as a preservative.

FOR *IN VITRO* LABORATORY USE ONLY.

**WARRANTY:** VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.