



**VMRD, Inc.**

PO Box 502, Pullman, WA 99163 USA

Telephone: + 1 (509) 334-5815

Fax: + 1 (509) 332-5356

E-mail: [vmrd@vmrd.com](mailto:vmrd@vmrd.com)

Web site: <http://www.vmrd.com>

## Certificate of Analysis

### **Porcine Circovirus (PCV)**

#### *Anti-Viral Antiserum*


**CATALOG NO.:** 210-70-PCRv

**VOLUME:** 2 ml

**LOT:** V01028-122006

**EXPIRATION:** 20 December 2006

**VIRUS:** Porcine Circovirus (PCV)



**DESCRIPTION:** Anti-PCV polyclonal antiserum. Liquid. Porcine origin. Made against a field isolate of type II; does not stain type I.

**INTENDED USE:** Useful for IFA.

**QUALITY CONTROL METHOD:** Indirect FA using VMRD PCV 10-well Substrate Slides (catalog no. 210-88-10-PCrV); Anti-Pig IgG FITC Conjugate (catalog no. 085-10).

**Specific Reaction:** 4+ positive fluorescence at 1:100.

**Other Reactions or Comments:** No background. Working dilution 1:800. IFA endpoint titer greater than 1:5,000.

**STORAGE:** This antiserum is provided in liquid form and should be stored at 4-8°C. DO NOT FREEZE! If antiserum becomes cloudy, it should be discarded. This antiserum contains 0.09% sodium azide as a preservative.

FOR IN VITRO LABORATORY USE ONLY.

**WARRANTY:** VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.