


Canine Distemper Virus (CDV)

Positive Blood Smear Slide

Catalog No.:	SLD-BSP-CDV
Shelf Life:	4 years from date of qualification
Substrate:	Blood positive for Canine Distemper Virus (CDV)
Species of Origin:	Canine



Description:

Canine blood smeared on glass slide positive for CDV antigen.

Intended Use:

As a fluorescent antigen positive control slide. Not intended for routine use as a control to show that the conjugate works. Rather, the CDV 2-well FA Control Slide (catalog no. SLD-FAC-CDV) should be used.

Storage:

Store sealed in foil pouch at -20°C. Avoid self-defrosting freezers.

References:

NA

Technical Data Sheet Version:

Version 1

FOR *IN VITRO* LABORATORY USE ONLY.

H:\Quality VMRD\COAs and TDS\TDS\Blood Smear\SLD-BSP-CDV TDS V1 CPM Active MKT.docx

23 September 2014

