



**VMRD**

PO Box 502, Pullman, WA 99163 USA

Telephone: 509-334-5815  
Fax: 509-332-5356  
E-mail: [vmrd@vmrd.com](mailto:vmrd@vmrd.com)  
Web site: <http://www.vmr.com>

---

## Certificate of Analysis

### **FA CONJUGATE DILUTING BUFFER**

**CATALOG NO.:** 210-91-CB  
**VOLUME:** 100 ml  
**LOT:** P080804-002  
**EXPIRATION:** 15 August 2010

**DESCRIPTION:** Conjugate Diluting Buffer is PBS with 1% BSA, preserved with 0.09% sodium azide.

**QUALITY CONTROL METHOD:** CDV, *Ehrlichia canis*, *Neospora caninum*, BVD, FIP-1, FIP-2, ERV/EHV, and CPV 12-well substrate slides tested with corresponding positive sera and anti-immunoglobulin FITC conjugates.

**Specific Reaction:** Positive sera dilute to same signal as previous lot.

**Other Comment:** Zero background.

**INTENDED USE:** Conjugate Diluting Buffer should be used for diluting concentrated conjugates for direct and indirect fluorescence.

**STORAGE:** Store at 4-8°C. DO NOT FREEZE!

FOR *IN VITRO* LABORATORY USE ONLY.

**WARRANTY:** VMRD, Inc. warrants that this product is as described in the quantity and contents stated on the label at the time of delivery to the customer. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE BEYOND THE LABEL DESCRIPTION, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. Remedy is limited to replacement of the product or refund of the purchase price. VMRD, Inc. is not liable for property damage, personal injury, or economic loss caused by the product. The information listed in this information sheet is provided for reference only, and should not be substituted for the user's own incoming material quality control.

